


EXHIBIT D

KS-DB60 R 32

clones 1, 2, 3, 4



Eco RI λ cut of insert



1	1	2	2	new	3	3	4	4
uncut	cut	uncut	cut	1ug	uncut	cut	uncut	cut
2ul	5ul	2ul	5ul		2ul	5ul	2ul	5ul

I picked this clone to grow up + subclone into H13

Southern:
Probed w/
 λ gt11 DB60
fragment



bel

Nitrocellulose
Blot


DB60 R32 - clone selected from Hare & DNH

λ gt 11 library of screening w/ DB60
EcoRI fragment - picked in 1st screen
w/ 2 P124F antibody

insert cut from λ gt 11, ~~for~~ subcloned
into KS⁻ bluescript plasmid

grown up in 50 ml of NM522 cells,
plasmid prep of 50 ml, run 10 μ g
on low salt gel + cut out EcoRI
fragment. Run over Clutip column
(very hard to push low salt wash through,
too much agarose, need smaller slice)

Fragment run on gel w/ marker



Fragment size is